

SOFTWARE PACKAGE
for
MEDICAL PHYSICS DEPARTMENT

This project is submitted in partial
fulfilment of the award
of the
DIPLOMA IN COMPUTER STUDIES
of the
HIGHER TECHNICAL INSTITUTE

CS/137

Project Supervisor : Mrs Eliza Angelidou
MSc BSc
HTI Lecturer, Computer Studies Course

Internal Accessor : Mr Stelios Christofides
PhD, CPhys, MInstP, MIPS M
Medical Physicist
Medical Physics Department
Nicosia General Hospital

External Accessor : Mr Yiannis Fouskotos
MA in Computer Science / Engineering

Design by

MARGARITA MAURAVOUNA

JUNE 1995




TABLE OF CONTENTS

1. INTRODUCTION.....	1
2. IMPORTANT INSTRUCTIONS	2
3. STARTING THE APPLICATION.....	3
3.1 INSTALLATION PROCEDURE.....	3
3.2 STARTING THE SYSTEM.....	4
4. PLOT GRAPH	7
4.1 SAVING A GRAPH	8
4.2 PRINTING THE CURVES	9
4.3 DELETING PEAKS	10
4.4 MODIFYING PEAKS	11
4.5 EXTRACTING GRAPHS.....	12
4.6 GENERATING RESULTS	13
4.7 ADDING GRIDLINES.....	14
4.8 REMOVING GRIDLINES	15
4.9 GETTING GRAPH INFORMATION	16
5. RETRIEVING A GRAPH	17
6. UTILITIES.....	18
6.1 EDITING BACKGROUND FILE	19
6.2 MAINTAINING USERS	20
6.2.1 ADDING A USER.....	21
6.2.2 DELETING A USER.....	22
6.2.3 MODIFYING A USER	23
6.3 BACKUP DATA FILES.....	24
6.4 RESTORE DATA FILES	25
7. GETTING GLOW CURVE INFORMATION.....	26

1. Introduction

The Glow Curve System handles data from TLD Readers 4000/6600 in a user friendly environment and also offers some useful utilities that are easy to use and save time as well.

Using the system is very simple even if you don't have any previous experience with computers.

You must refer to this manual every time you have some difficulties .

Everything about the operations of the system is explained in this user guide.